



**BIOGENETS**

## NOMBRE DEL PRODUCTO

### Analizador de Inmunoensayo Fluorescente – Modelos BJC-FT

Ofrecemos una amplia gama de productos por inmunofluorescencia.

#### Analyzer

Product	Sample	Cut-Off	Test Range	Pack
Fluorescence Immunoassay Analyzer version 1.0 (without scanner)	/	/	/	1pc
Fluorescence Immunoassay Analyzer version 1.0 (with scanner)	/	/	/	1pc
Incubox for HBA1C Test (Temperature incubator)	/	/	/	1pc

#### Inflamation Marker

Product	Sample	Cut-Off	Test Range	Pack
CRP Test Cassette	WB/S/P	1.0 ug/mL, 10ug/mL	0.5~200 mg/L	10T/25T
CPT Test Cassette	WB/S/P	0.5 ng/mL	0.1~50 ng/ml	10T/25T
FOB Test Cassette	Feces	50ng/ml	10~ 1000ng/ml	10T/25T

#### Tumor Marker

Product	Sample	Cut-Off	Test Range	Pack
PSA Test Cassette	S/P	4 ng/mL	2~100 ng/mL	10T/25T
AFP Test Cassette	S/P	20 ng/mL	2~100 ng/mL	10T/25T
CEA Test Cassette	S/P	4ng/ml	1~500 ng/mL	10T/25T



**BIOGENETS**

### Cardiac Marker

Product	Sample	Cut-Off	Test Range	Pack
cTnI Test Cassette	W/S/P	0.5 ng/mL	0.2-40 ng/mL	10T/25T
cTnT Test Cassette	W/S/P	0.5ng/ml	0.2-40ng/ml	10T/25T
NT-proBNP Test Cassette	W/S/P	< 300 pg/mL	18-35000 pg/mL	10T/25T
CK-MB Test Cassette	W/S/P	< 5 ng/mL	0.3-100 ng/mL	10T/25T
H-FABP	W/S/P	< 8ng/mL	1-120 ng/mL	10T/25T
Myoglobin Test Cassette	W/S/P	< 90 ng/mL	5-200 ng/mL	10T/25T

### Coagulation

Product	Sample	Cut-Off	Test Range	Pack
D-Dimer Test Cassette	W/S/P	0.5 mg/L	0.1~10 mg/L	10T/25T

### Endocrinology

Product	Sample	Cut-Off	Test Range	Pack
FSH Test Cassette	W/S/P	≤ 25mIU/mL	10~80mIU/mL	10T/25T
LH Test Cassette	W/S/P	≤ 40mIU/mL	20 ~300mIU/mL	10T/25T
T4 Test Cassette	S/P	65~155nmol/L	12.87 ~300nmol	10T/25T
β-HCG Test Cassette	S/P	< 25mIU/ml	5~10,000mIU /m	10T/25T
TSH Test Cassette	S/P	< 5 μIU/mL	0.1- 100μIU/mL	10T/25T
T3 Test Cassette	S/P	0.8-	0.62 ~9.24nmol/	10T/25T
Testosterone Test Cassette	S/P	/	1-20 ng/mL	10T/25T
Progesterone(P4) Test Cassette	S/P	/	1.4-60 ng/mL	10T/25T



**BIOGENETS**

### Nephrology

Product	Sample	Cut-Off	Test Range	Pack
mAlb Test Cassette	Urine	< 20mg/L	5-300mg/L	10T/25T
CysC Test Cassette	W/S/P	0.51-1.09mg/L	0.5-10mg/L	10T/25T
β2MG Test Cassette	W/S/P	≤ 2.2mg/l	0.2~20mg/L	10T/25T
N-GAL Test Cassette	Urine	< 13.6 ng/mL	10-1500ng/ml	10T/25T

### Diabetes

Product	Sample	Cut-Off	Test Range	Pack
HbA1c Test Cassette	WB	4-6%	4-14.5%	10T/25T

### Infectious Diseases

Product	Sample	Cut-Off	Test Range	Pack
Strep A Test Cassette	Swab	/	Qualitative	10T/20T
FLU A+B Test Cassette	Swab	/	Qualitative	10T/20T
RSV Test Cassette	W/S/P	/	Qualitative	10T/20T
Dengue IgG/IgM Test Cassette	W/S/P	/	Qualitative	10T/20T
Dengue NS1 Test Cassette	W/S/P	/	Qualitative	10T/20T
COVID-19 antibody Test Cassette	W/S/P	/	Qualitative	10T/20T
COVID-19 antigen Test Cassette	Swab	/	Qualitative	10T/20T

### Autoimmunity

Product	Sample	Cut-Off	Test Range	Pack
IgE Test Cassette	W/S/P	330IU/ml	20 ~1600IU/mL	10T/25T



**BIOGENETS**

## NOMBRE DEL PRODUCTO

### Analizador de Inmunoensayo Fluorescente – Modelos BJC-NT

Ofrecemos una amplia gama de productos por inmunofluorescencia.

#### Analyzer

Product	Sample	Cut-Off	Test Range	Pack
AFR-200 Fluorescence Immunoassay Analyzer 2.0 (without scanner) / / / 1 pc CE	/	/	/	1pc
Fluorescence Immunoassay Analyzer version 2.0 (with scanner)	/	/	/	1pc
Incubox for HBA1C Test (Temperature incubator)	/	/	/	1pc

#### Inflammation Marker

Product	Sample	Cut-Off	Test Range	Pack
CRP Test Cassette	WB/S/P	1.0 ug/mL, 10ug/mL	0.5~200 mg/L	10T/25T
PCT Test Cassette	WB/S/P	0.5 ng/mL	0.1~50 ng/ml	10T/25T
FOB Test Cassette	Feces	50ng/ml	10~ 1000ng/ml	10T/25T

#### Tumor Marker

Product	Sample	Cut-Off	Test Range	Pack
PSA Test Cassette	S/P	4 ng/mL	2~100 ng/mL	10T/25T
AFP Test Cassette	S/P	20 ng/mL	5~400 ng/mL	10T/25T
CEA Test Cassette	S/P	4ng/ml	1~500 ng/mL	10T/25T



**BIOGENETS**

**Cardiac Marker**

<b>Product</b>	<b>Sample</b>	<b>Cut-Off</b>	<b>Test Range</b>	<b>Pack</b>
<b>cTnI Test Cassette</b>	S/P	0.5 ng/mL	0.1-40 ng/mL	<b>10T/25T</b>
<b>cTnT Test Cassette</b>	W/S/P	0.5ng/ml	0.2-40ng/ml	<b>10T/25T</b>
<b>NT-proBNP Test Cassette</b>	W/S/P	< 300 pg/mL	18-35000 pg/mL	<b>10T/25T</b>
<b>CK-MB Test Cassette</b>	W/S/P	< 5 ng/mL	0.3-100 ng/mL	<b>10T/25T</b>
<b>H-FABP</b>	W/S/P	< 8ng/mL	1-120 ng/mL	<b>10T/25T</b>
<b>Myoglobin Test Cassette</b>	<b>W/S/P</b>	<b>&lt; 90 ng/mL</b>	<b>5-200 ng/mL</b>	<b>10T/25T</b>

**Coagulation**

<b>Product</b>	<b>Sample</b>	<b>Cut-Off</b>	<b>Test Range</b>	<b>Pack</b>
<b>D-Dimer Test Cassette</b>	<b>W/S/P</b>	<b>0.5 mg/L</b>	<b>0.1~10 mg/L</b>	<b>10T/25T</b>

**Endocrinology**

<b>Product</b>	<b>Sample</b>	<b>Cut-Off</b>	<b>Test Range</b>	<b>Pack</b>
<b>FSH Test Cassette</b>	W/S/P	≤ 25mIU/mL	10~80mIU/mL	<b>10T/25T</b>
<b>LH Test Cassette</b>	W/S/P	≤ 40mIU/mL	20 ~300mIU/mL	<b>10T/25T</b>
<b>T4 Test Cassette</b>	S/P	65~155nmol/L	12.87 ~300nmol	<b>10T/25T</b>
<b>β-HCG Test Cassette</b>	S/P	< 25mIU/ml	5~10,000mIU /m	<b>10T/25T</b>
<b>TSH Test Cassette</b>	S/P	< 5 μIU/mL	0.1- 100μIU/mL	<b>10T/25T</b>
<b>T3 Test Cassette</b>	S/P	0.8- 2.0ng/mL (1.23 ~3.08nmol/L)	0.62 ~9.24nmol/	<b>10T/25T</b>
<b>Testosterone Test Cassette</b>	S/P	/	1-20 ng/mL	<b>10T/25T</b>
<b>Progesterone(P4) Test Cassette</b>	<b>S/P</b>	<b>/</b>	<b>1.4-60 ng/mL</b>	<b>10T/25T</b>



**BIOGENETS**

### Nephrology

Product	Sample	Cut-Off	Test Range	Pack
mAlb Test Cassette	Urine	< 20mg/L	5-300mg/L	10T/25T
CysC Test Cassette	W/S/P	0.51-1.09mg/L	0.5-10mg/L	10T/25T
β2MG Test Cassette	W/S/P	≤ 2.2mg/l	0.2~20mg/L	10T/25T
N-GAL Test Cassette	Urine	< 13.6 ng/mL	10-1500ng/ml	10T/25T

### Diabetes

Product	Sample	Cut-Off	Test Range	Pack
HbA1c Test Cassette	WB	4-6%	4-14.5%	10T/25T

### Infectious Diseases

Product	Sample	Cut-Off	Test Range	Pack
Strep A Test Cassette	Swab	Negative	≥1E6	10T/25T
FLU A+B Test Cassette	Swab	Negative	30HA/ml	10T/25T
RSV Test Cassette	Swab	/	Qualitative	10T/20T
Dengue IgG/IgM Test Cassette	W/S/P	/	Qualitative	10T/25T
Dengue NS1 Test Cassette	W/S/P	/	Qualitative	10T/25T
COVID-19 antibody Test Cassette	W/S/P	/	Qualitative	10T/25T
COVID-19 antigen Test Cassette	Swab	/	Qualitative	10T/20T

### Autoimmunity

Product	Sample	Cut-Off	Test Range	Pack
IgE Test Cassette	W/S/P	330IU/ml	20 ~1600IU/mL	10T/25T



**BIOGENETS**

## NOMBRE DEL PRODUCTO

### Analizador de Inmunoensayo Fluorescente – Modelos BJC-FTG

Ofrecemos una amplia gama de productos por inmunofluorescencia.

#### Analyzer

Product	Sample	Cut-Off	Test Range	Pack
Fluorescence Immunoassay Analyzer (Portable)	/	/	/	1pc
Incubox for HBA1C Test (Temperature incubator)	/	/	/	1pc

#### Inflammation Marker

Product	Sample	Cut-Off	Test Range	Pack
CRP Test Cassette	WB/S/P	1.0 ug/mL, 10ug/mL	0.5~200 mg/L	10T/25T
PCT Test Cassette	WB/S/P	0.5 ng/mL	0.1~50 ng/ml	10T/25T
FOB Test Cassette	Feces	50ng/ml	10~ 1000ng/ml	10T/25T

#### Tumor Marker

Product	Sample	Cut-Off	Test Range	Pack
PSA Test Cassette	S/P	4 ng/mL	2~100 ng/mL	10T/25T
AFP Test Cassette	S/P	20 ng/mL	5~400 ng/mL	10T/25T
CEA Test Cassette	S/P	4ng/ml	1~500 ng/mL	10T/25T



**BIOGENETS**

**Cardiac Marker**

<b>Product</b>	<b>Sample</b>	<b>Cut-Off</b>	<b>Test Range</b>	<b>Pack</b>
<b>cTnI Test Cassette</b>	S/P	0.5 ng/mL	0.1-40 ng/mL	<b>10T/25T</b>
<b>cTnT Test Cassette</b>	W/S/P	0.5ng/ml	0.2-40ng/ml	<b>10T/25T</b>
<b>NT-proBNP Test Cassette</b>	W/S/P	< 300 pg/mL	18-35000 pg/mL	<b>10T/25T</b>
<b>CK-MB Test Cassette</b>	W/S/P	< 5 ng/mL	0.3-100 ng/mL	<b>10T/25T</b>
<b>H-FABP</b>	W/S/P	< 8ng/mL	1-120 ng/mL	<b>10T/25T</b>
<b>Myoglobin Test Cassette</b>	<b>W/S/P</b>	<b>&lt; 90 ng/mL</b>	<b>5-200 ng/mL</b>	<b>10T/25T</b>

**Coagulation**

<b>Product</b>	<b>Sample</b>	<b>Cut-Off</b>	<b>Test Range</b>	<b>Pack</b>
<b>D-Dimer Test Cassette</b>	<b>W/S/P</b>	<b>0.5 mg/L</b>	<b>0.1~10 mg/L</b>	<b>10T/25T</b>

**Endocrinology**

<b>Product</b>	<b>Sample</b>	<b>Cut-Off</b>	<b>Test Range</b>	<b>Pack</b>
<b>FSH Test Cassette</b>	W/S/P	≤ 25mIU/mL	10~80mIU/mL	<b>10T/25T</b>
<b>LH Test Cassette</b>	W/S/P	≤ 40mIU/mL	20 ~300mIU/mL	<b>10T/25T</b>
<b>T4 Test Cassette</b>	S/P	65~155nmol/L	12.87 ~300nmol	<b>10T/25T</b>
<b>β-HCG Test Cassette</b>	S/P	< 25mIU/ml	5~10,000mIU /m	<b>10T/25T</b>
<b>TSH Test Cassette</b>	S/P	< 5 μIU/mL	0.1- 100μIU/mL	<b>10T/25T</b>
<b>T3 Test Cassette</b>	S/P	0.8- 2.0ng/mL (1.23 ~3.08nmol/L)	0.62 ~9.24nmol/	<b>10T/25T</b>
<b>Testosterone Test Cassette</b>	S/P	/	1-20 ng/mL	<b>10T/25T</b>
<b>Progesterone(P4) Test Cassette</b>	<b>S/P</b>	<b>/</b>	<b>1.4-60 ng/mL</b>	<b>10T/25T</b>





**BIOGENETS**

### Nephrology

Product	Sample	Cut-Off	Test Range	Pack
mAlb Test Cassette	Urine	< 20mg/L	5-300mg/L	10T/25T
CysC Test Cassette	W/S/P	0.51-1.09mg/L	0.5-10mg/L	10T/25T
β2MG Test Cassette	W/S/P	≤ 2.2mg/l	0.2~20mg/L	10T/25T
N-GAL Test Cassette	Urine	< 13.6 ng/mL	10-1500ng/ml	10T/25T

### Diabetes

Product	Sample	Cut-Off	Test Range	Pack
HbA1c Test Cassette	WB	4-6%	4-14.5%	10T/25T

### Infectious Diseases

Product	Sample	Cut-Off	Test Range	Pack
Strep A Test Cassette	Swab	Negative	≥1E6	10T/25T
FLU A+B Test Cassette	Swab	Negative	30HA/ml	10T/25T
RSV Test Cassette	Swab	/	Qualitative	10T/20T
Dengue IgG/IgM Test Cassette	W/S/P	/	Qualitative	10T/25T
Dengue NS1 Test Cassette	W/S/P	/	Qualitative	10T/25T
COVID-19 antibody Test Cassette	W/S/P	/	Qualitative	10T/25T
COVID-19 antigen Test Cassette	Swab	/	Qualitative	10T/20T

### Autoimmunity

Product	Sample	Cut-Off	Test Range	Pack
IgE Test Cassette	W/S/P	330IU/ml	20 ~1600IU/mL	10T/25T