



BIOGENETS

NOMBRE DEL PRODUCTO

Marcadores Cardiacos

A continuación, se encuentra listado nuestro panel completo de marcadores cardiacos en formato de test rápidos individuales y en combo.

Product Description	Specimen	Format	Kit Size	Cut-Off
Cardiac Troponin I (cTnI) Rapid Test Cassette	WB/S/P	Cassette	25 T	0.5 ng/mL
Cardiac Troponin T (cTnT) Rapid Test Cassette	WB/S/P	Cassette	10 T	0.5 ng/mL
CK-MB Rapid Test Cassette	WB/S/P	Cassette	25 T	5 ng/mL
Myoglobin Rapid Test Cassette	WB/S/P	Cassette	25 T	50 ng/mL
H-FABP Rapid Test Cassette	WB/S/P	Cassette	10 T	8 ng/mL
H-FABP and cTnI Combo Rapid Test Cassette	WB/S/P	Cassette	10 T	See Insert
Myoglobin/CK-MB/Troponin I Combo Rapid Test Cassette	WB/S/P	Cassette	25 T	See Insert
H-FABP and Myoglobin/CK-MB/Cardiac Troponin I Combo Rapid Test Cassette	WB/S/P	Cassette	10 T	See Insert
NT-proBNP Rapid Test Cassette	WB/S/P	Cassette	10 T	0.45 ng/mL
D-dimer Rapid Test Cassette	WB / P	Cassette	10 T	500 ng/mL
C-reactive protein Rapid Test Dipstick	WB/S/P	Dipstick	20 T	10 µg/mL
C-reactive protein Rapid Test Cassette	WB/S/P	Cassette	20 T	10 µg/mL*, 0.5 µg/mL
C-reactive protein Semi-Quantitative Rapid Test Dipstick	WB/S/P	Dipstick	20 T	10-40-80 µg/mL

C-reactive protein Semi-Quantitative Rapid Test Cassette	WB/S/P	Cassette	20 T	10-30 µg/mL
C-reactive protein Semi-Quantitative Rapid Test Cassette	WB/S/P	Cassette	20 T	10-40-80 µg/mL
Procalcitonin(PCT) Rapid Test Cassette	S / P	Cassette	10 T	1 ng/mL
Procalcitonin(PCT) Rapid Test Cassette	WB/S/P	Cassette	10 T	1 ng/mL