



B I O G E N E T S

NOMBRE DEL PRODUCTO

Test Rápidos de Oncología

A continuación, se encuentra listado nuestro panel completo de marcadores tumorales en formato de test rápidos.

Product Description	Specimen	Format	Kit Size	Cut-Off
AFP Rapid Test Cassette	WB/S/P	Cassette	25 T	10 ng/mL
CEA Rapid Test Cassette	WB/S/P	Cassette	25 T	5 ng/mL
PSA Rapid Test Cassette	WB/S/P	Cassette	25 T	3-10 ng/mL
PSA Qualitative Rapid Test Dipstick	WB/S/P	Dipstick	25 T	4 ng/mL
PSA Qualitative Rapid Test Cassette	WB/S/P	Cassette	25 T	4 ng/mL
CA125 Rapid Test Cassette	WB/S/P	Cassette	10 T	40 U/mL
CA15-3 Rapid Test Cassette	WB/S/P	Cassette	10 T	30 U/mL
CA19-9 Rapid Test Cassette	WB/S/P	Cassette	10 T	40 U/mL
FOB Rapid Test Cassette	Feces	Cassette	25 T	40 ng/mL, 50 ng/mL*
FOB Rapid Test Cassette	Feces	Cassette	1 T	40 ng/mL, 50 ng/mL*
FOB Rapid Test Dipstick	Feces	Dipstick	25 T	40 ng/mL, 50 ng/mL*
FOB Rapid Test Dipstick	Feces	Dipstick	1 T	40 ng/mL, 50 ng/mL*
Calprotectin and FOB Combo Rapid Test Cassette	Feces	Cassette	10 T	See Insert

Transferrin and FOB Combo Rapid Test Cassette	Feces	Cassette	25 T	See Insert
Hb+Hb-Hp Combo Rapid Test Cassette	Feces	Cassette	25 T	50 ng/mL
Transferrin/FOB and Hb-Hp Combo Rapid Test Cassette	Feces	Cassette	25 T	See Insert